



## Neuroform Atlas™




### STENT SYSTEM

### Device Comparison Guide

with LVIS Jr. and Leo Baby

#### Purpose

To present technical data obtained in bench tests involving three stents approved for treating wide-necked aneurysms.

	Neuroform Atlas Stent	MicroVention LVIS Jr. Stent	Balt Extrusion Leo Baby Stent
<b>Design</b>	Open Cell	Woven (12-wire)	Woven (16-wire)
<b>Recommended vessel diameter (mm)</b>	2.0-4.5	2.0-3.0	1.5-3.1
<b>Available Stent Diameter(s) (mm)</b>	3.0, 4.0, 4.5	3.5	2.0, 2.5
<b>Available Stent Lengths (mm)</b>	15, 21, 24, 30	18, 23, 28, 33	12, 18, 25
<b>Minimum Microcatheter ID (in)</b>	0.0165	0.0170	0.0165
<b>Strut Width</b>	57 microns (8-crown) 31 microns (12-crown)	54 micron wire	40 micron wire
<b>Strut Thickness</b>	55 microns	54 micron wire	40 micron wire
<b>Radiopacity Markers</b>	3 Platinum/Iridium Alloy markers at each end	3 Tantalum markers at each end and 3 longitudinal helical Tantalum strands	2 longitudinal helical Platinum coils
<b>Metal-to-Artery Ratio</b>	6-12%	12-16%	12-14%
<b>Foreshortening</b>	2.9-6.3%	8.0-40.0%	16.1-32.3%
<b>Stent Delivery Wire Tip Length (mm)</b>	0.0	5.9	8.0
<b>Device Image</b>			

Bench test results. n=1. Bench test results may not necessarily be indicative of clinical performance. Testing performed by Stryker Neurovascular. Data are on file at Stryker Neurovascular and will be made available upon request. All photographs are taken by Stryker Neurovascular.

Information for LVIS Jr. excerpted from the literature published by MicroVention, LVIS Jr. IFU - PD070019D Revised 2012-05.

Information for Leo Baby excerpted from the literature published by Balt Extrusion, DCO15GB Revised 2012/05.

**Neuroform Atlas™  
Stent System**

**See package insert for complete indications, contraindications, warnings and instructions for use.**

**INTENDED USE/INDICATIONS FOR USE**

The Neuroform Atlas Stent System is intended to be used with occlusive devices in the treatment of intracranial aneurysms.

**THIS DOCUMENT IS INTENDED SOLELY FOR THE USE OF HEALTHCARE PROFESSIONALS.**

A physician must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that physicians be trained in the use of any particular product before using it in a procedure. The information presented is intended to demonstrate the breadth of Stryker product offerings. A physician must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. The Stryker products listed above are CE marked according to the Medical Device Directive 93/42/EEC.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Neuroform Atlas. All other trademarks are trademarks of their respective owners or holders.

Copyright © 2015 Stryker  
NV00016425.AA | Page 2 of 2



RAQA Manager  
**Stryker France S.A.S.**  
ZAC-Avenue de Satolas Green  
69330 Pusignan  
France



**Stryker Neurovascular**  
47900 Bayside Parkway  
Fremont, CA 94538

[strykerneurovascular.com](http://strykerneurovascular.com)

Date of Release: AUG/2015

EX\_EN\_IL